

ALMORE

microfil composite instruments

Lightweight composite instruments are ideal for composite restoration development. The instruments are autoclavable and the tips at end of the instruments (with black plastic) are replaceable. The instruments unique thin and flexible ends are particularly helpful in shaping light-cured composites. The thin blades easily slide under the gum line to shape the margin. The flat surfaces contour flat facial surfaces of the teeth and the thin edge is used for adding esthetic highlights to the tooth's surface. **DO NOT PLACE IN ULTRASONIC!**



96041 GOLD INSTRUMENT: The replaceable end's spatula shape is used on the facial surfaces of the maxillary incisors and the buccals of the posterior teeth or where the surface is large and rounded. It is excellent for adapting a gum line margin, and for forming the gentle curve toward the interproximal and the edge for esthetic highlights.

The gold instrument's permanent end has a sturdy curved spoon shape. The back of the instrument is used for shaping the curves of the lingual concavities. The front side of the spoon is used to form the curvature at the gum line or interproximal on the lingual side of the teeth and also may be used to carry the material to the preparation.



96043 BLUE INSTRUMENT: The replaceable end's spear shape cutting edge on the blue instrument is used for trimming the composite material especially at the margin under the gum or interproximals.

The blue instrument's permanent end is especially designed to reach between the teeth under the gums. The small pointed tip is perfect for adapting the gingival margin interproximally. It also has a cutting edge on the inside of the sickle which can be used to trim interproximally or to shave the facial surface of the tooth if it is laid flat against the surface and gently moved across the tooth plane. It is also helpful for shaping the proximal lingual angles of the teeth when it is placed against the surface, using the back of the sickle and moving it toward the interproximal of the tooth.



96042 GREEN INSTRUMENT: The spear-shape replaceable end is used in the manner as the Gold instrument's spatula, but is used for smaller teeth such as the lower anterior teeth and is especially helpful for creating the esthetic highlights of the proximal line angles.

The permanent end of the Green instrument has a deeper more flexible spoon shape and is used where more contour is needed.

96044 SET OF ALL (3) INSTRUMENTS (GOLD, GREEN, BLUE)